

*Experience the Difference*

## ISBT Blood Bank Labels



**TimeMed Labeling Systems, Inc.**

144 Tower Drive  
Burr Ridge, IL 60527  
Phone: 800-323-4840  
Fax: 800-548-5359  
[www.TimeMed.com](http://www.TimeMed.com)

**PHARMEX**<sup>®</sup>  
A Division of TimeMed

**T**imeMed can help you meet the ISBT Blood Bank standards head-on with our complete line of blank and pre-printed labels.

- ♦ Precisely sized for on-demand printing
- ♦ Quality labels that easily run through your thermal printer
- ♦ Produce accurate bar codes and eye readable information
- ♦ Meets U.S. Food and Drug Administration requirements

Contact your TimeMed Sales Representative at 800-323-4840 for additional information on labels and thermal transfer ribbons that are compatible with your ISBT 128 on-demand printing systems.

# ISBT Blood Bank Labels

Our on-demand, blank stock labels are available in thermal transfer and direct thermal materials to help you meet the specific needs of your facility's blood bank.

## Blank Stock ISBT Labels

Product Code	Label Size	Core Size	Material/Adhesive	Qty.
TD-128-44	4" x 4"	3"	Coated Thermal	950/RL - 4RLS/BX
TDS-128-44	4" x 4"	3"	Coated Thermal Synthetic	950/RL - 4RLS/BX
TTR-128-44	4" x 4"	3"	Thermal Transfer	950/RL - 4RLS/BX
TTRS-128-44	4" x 4"	3"	Thermal Transfer Synthetic	950/RL - 4RLS/BX

## Stock Thermal Transfer Ribbons

Product Code	For Use With	Width Millimeter	Width Inches	Length Meter	Length Feet	Box Qty.
<b>Zebra Z4M</b>						
TTR-FRD10275	TTR	102	4.02	450	1476	6
TTR-FRD10238	TTR/TTRS	102	4.02	450	1476	6
TTR-FRD10233	TTRS	102	4.02	450	1476	1
<b>Sato BB-412e</b>						
TTR-FST10275	TTR	102	4.02	450	1345	6
TTR-FST10208	TTR/TTRS	102	4.02	450	1476	6
TTR-FST10230	TTRS	102	4.02	450	1476	6

If you're looking to be compliant with ISBT 128, let TimeMed assist you with preprinted labels like sequential numbers for donor identification, ABO/Rh labels, and more. Visit [www.TimeMed.com](http://www.TimeMed.com) to view the most recent additions to our preprinted ISBT label product line.

## Pre-printed ABO/Rh Labels

Product Code	Desc.	Label Size	Core Size	Material/Adhesive	Qty.
ISBT0600	A Rh -	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX
ISBT6200	A Rh +	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX
ISBT1700	B Rh -	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX
ISBT7300	B Rh +	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX
ISBT9500	O Rh -	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX
ISBT5100	O Rh +	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX
ISBT2800	AB Rh -	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX
ISBT8400	AB Rh +	2" x 2"	3"	Thermal Transfer Synthetic	250/RL - 1RL/BX

