

Experience the Difference

StainerShield® Slide Identification Labels



TimeMed®

TimeMed Labeling Systems, Inc.

144 Tower Drive
Burr Ridge, IL 60527
Phone: 800-323-4840
Fax: 800-548-5359
www.TimeMed.com

PHARMEX®
A Division of TimeMed





Tired of spending valuable time and money on misidentified specimens, redraw samples and repeat testing?

TimeMed has an effective solution for reliable pre-stainer slide labeling. StainerShield® slide labels are a simplified solution to the complex problem of properly identifying and tracking specimen slides in a laboratory environment. Lab technicians can now print on-demand, patient-specific labels and apply them to tissue and specimen slides as the slides are produced, before they go through the stainer. No overlamine flap, UV ink, or etching system is needed.

Are you looking for a way to error proof your process and streamline your workflow to eliminate costly inefficiencies?


These labels can be printed on-demand using most standard desktop thermal printers. In one simple printing step, technicians can print slide labels with both text and bar codes – including 2D bar codes – that contain data directly from the lab's information system.

Visit our website at www.TimeMed.com for more information on our reliable pre-stainer labeling solution. *See the reverse side for a variety of sizes, colors and configurations of StainerShield® slide labels.*






1-Up Style (Click on a Product Code link below to place an order)

Product Code	Label Size	Core Size	Max OD	Mat/Adh	Qty.
TDSS1-1-1278	1/2" x 7/8"	1"	4-1/5"	DT	1.8M/RL-4RLS/BX
TDSS1-1-7878	7/8" x 7/8"	1"	4-13/32"	DT	1.8M/RL-4RLS/BX
TDSS15-17878	7/8" x 7/8"	1-1/2"	4-13/32"	DT	1.8M/RL-4RLS/BX
TDSS3-1-7878	7/8" x 7/8"	3"	5"	DT	1.8M/RL-4RLS/BX
TDSS1-17878P	7/8" x 7/8"	1"	4-13/32"	DT	1.8M/RL-4RLS/BX
TDSS1-17878B	7/8" x 7/8"	1"	4-13/32"	DT	1.8M/RL-4RLS/BX
TDSS1-17878Y	7/8" x 7/8"	1"	4-13/32"	DT	1.8M/RL-4RLS/BX




2 Across Style (Click on a Product Code link below to place an order)

Product Code	Label Size	Core Size	Max OD	Mat/Adh	Qty.
TDSS1-2-7878	7/8" x 7/8"	1"	4-1/5"	DT	3.6M/RL-2RLS/BX
TDSS3-2-7878	7/8" x 7/8"	3"	5"	DT	3.6M/RL-2RLS/BX



4 Across Style (Click on a Product Code link below to place an order)

Product Code	Label Size	Core Size	Max OD	Mat/Adh	Qty.
TDSS1-4-7878	7/8" x 7/8"	1"	4-1/5"	DT	7.2M/RL-1RL/BX
TDSS1-47878B	7/8" x 7/8"	1"	4-1/5"	DT	7.2M/RL-1RL/BX
TDSS3-4-7878	7/8" x 7/8"	3"	5"	DT	7.2M/RL-1RL/BX
TDSS1-4-1278	1/2" x 7/8"	1"	4-1/5"	DT	7.2M/RL-1RL/BX
TDSS3-4-1278	1/2" x 7/8"	3"	4-1/5"	DT	7.2M/RL-1RL/BX



For use with Cerner Millennium® Systems (Click on a Product Code link below to place order)

Product Code	Label Size	Core Size	Max OD	Mat/Adh	Qty.
TDSS1-4-1516	15/16" x 15/16"	1"	4-1/5"	DT	7.2M/RL-1RL/BX
TDSS3-4-1516	15/16" x 15/16"	3"	4-1/5"	DT	7.2M/RL-1RL/BX

Not Just Another "Xylene-resistant" Label!

The labels' construction utilizes unique and patented polymeric technology that is specifically designed to withstand laboratory specimen processing. Once applied to the slide, StainerShield® slide labels are designed to resist a stainer's harsh reagents and solvents without fading, degrading or separating from the slide. Both the text and the bar codes printed on the label remain legible and produce successful scans, even after stainer processing.

Additional StainerShield® labels are available in: yellow, pink & light blue. Contact your TimeMed Sales Representative for more details.



"Slides must be adequately identified. . .adequate procedures must be in place to ensure positive identification of slides during processing. The permanent slide label must be legible and indelible.

- CAP Laboratory Accreditation Checklist ANP21150